

## Product datasheet for **TA341240**


### ATP1B3 Rabbit Polyclonal Antibody

#### Product data:

Product Type:	Primary Antibodies
Applications:	IHC
Recommended Dilution:	IHC-P (5 µg/ml)
Reactivity:	Human
Host:	Rabbit
Clonality:	Polyclonal
Immunogen:	ATP1B3 antibody was raised against synthetic 18 amino acid peptide from internal region of human ATP1B3. Percent identity with other species by BLAST analysis: Human, Monkey (100%); Gibbon (94%).
Formulation:	PBS, 0.1% sodium azide.
Concentration:	lot specific
Purification:	Immunoaffinity purified
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Gene Name:	ATPase Na <sup>+</sup> /K <sup>+</sup> transporting subunit beta 3
Database Link:	<a href="#">NP_001670</a> <a href="#">Entrez Gene 483 Human</a> <a href="#">P54709</a>
Synonyms:	ATPB-3; CD298
Protein Families:	Transmembrane
Protein Pathways:	Cardiac muscle contraction



[View online »](#)

**Product images:**

Anti-ATP1B3 antibody IHC of human prostate. Immunohistochemistry of formalin-fixed, paraffin-embedded tissue after heat-induced antigen retrieval.