

Product datasheet for **TA341243**


ARTS1 (ERAP1) Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	IHC
Recommended Dilution:	IHC-P (2.5 µg/ml)
Reactivity:	Human (Predicted: Monkey, Bat, Dog, Pig, Rabbit)
Host:	Rabbit
Clonality:	Polyclonal
Immunogen:	ARTS1 / ERAP1 antibody was raised against synthetic 15 amino acid peptide from C-Terminus of human ERAP1 / ARTS1. Percent identity with other species by BLAST analysis: Human, Chimpanzee, Gorilla (100%); Monkey, Galago, Marmoset, Elephant, Dog, Bat, Rabbit, Pig (93%); Orangutan, Gibbon, Mouse, Rat, Panda, Bovine, Horse, Guinea pig (87%).
Formulation:	PBS, 0.1% sodium azide.
Concentration:	lot specific
Purification:	Immunoaffinity purified
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Gene Name:	endoplasmic reticulum aminopeptidase 1
Database Link:	NP_057526 Entrez Gene 488898 Dog Entrez Gene 699996 Monkey Entrez Gene 51752 Human Q9NZ08
Synonyms:	A-LAP; ALAP; APPILS; ARTS-1; ARTS1; ERAAP; ERAAP1; PILS-AP; PILSAP
Protein Families:	Druggable Genome, Protease, Secreted Protein



[View online »](#)

Product images:

Anti-ERAP1 / ARTS1 antibody IHC of human kidney. Immunohistochemistry of formalin-fixed, paraffin-embedded tissue after heat-induced antigen retrieval.