

## Product datasheet for **TA341311**


### GM2A Rabbit Polyclonal Antibody

#### Product data:


Product Type:	Primary Antibodies
Applications:	IHC
Recommended Dilution:	IHC-P (10 - 12.5 µg/ml)
Reactivity:	Human
Host:	Rabbit
Clonality:	Polyclonal
Immunogen:	GM2A antibody was raised against synthetic 14 amino acid peptide from N-terminus of human GM2A. Percent identity with other species by BLAST analysis: Human, Gorilla, Gibbon (100%); Monkey (86%).
Formulation:	PBS, 0.1% sodium azide.
Concentration:	lot specific
Purification:	Immunoaffinity purified
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Gene Name:	GM2 ganglioside activator
Database Link:	<a href="#">NP_000396</a> <a href="#">Entrez Gene 2760 Human</a> <a href="#">P17900</a>
Synonyms:	GM2-AP; SAP-3
Protein Families:	Druggable Genome
Protein Pathways:	Lysosome



[View online »](#)

**Product images:**

Anti-GM2A antibody IHC of human kidney.  
Immunohistochemistry of formalin-fixed,  
paraffin-embedded tissue after heat-induced  
antigen retrieval.



Anti-GM2A antibody IHC of human skin.  
Immunohistochemistry of formalin-fixed,  
paraffin-embedded tissue after heat-induced  
antigen retrieval.