


Product datasheet for **TA341473**

ZNF641 Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	WB
Recommended Dilution:	WB
Reactivity:	Human
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	The immunogen for anti-ZNF641 antibody: synthetic peptide directed towards the middle region of human ZNF641. Synthetic peptide located within the following region: CGIVSLRFPIPKLDMLSQLEGGEQWVPDPQDLEERDILRVTYTGGGSE
Formulation:	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Concentration:	lot specific
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Predicted Protein Size:	49 kDa
Gene Name:	zinc finger protein 641
Database Link:	NP_689533 Entrez Gene 121274 Human Q96N77
Background:	ZNF641 is a transcriptional activator. It activates transcriptional activities of SRE and AP-1.
Synonyms:	DKFZp667D1012; FLJ31295
Note:	Immunogen Sequence Homology: Dog: 100%; Pig: 100%; Rat: 100%; Horse: 100%; Human: 100%; Rabbit: 100%; Guinea pig: 93%; Bovine: 92%; Mouse: 79%
Protein Families:	Transcription Factors


[View online »](#)

Product images:


WB Suggested Anti-ZNF641 Antibody Titration:
 0.2-1 ug/ml; ELISA Titer: 1:312500; Positive
 Control: 721_B cell lysate