

## Product datasheet for **TA341477**

### ZNF440 Rabbit Polyclonal Antibody

#### Product data:

Product Type:	Primary Antibodies
Applications:	WB
Recommended Dilution:	WB
Reactivity:	Human
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	The immunogen for anti-ZNF440 antibody: synthetic peptide directed towards the C terminal of human ZNF440. Synthetic peptide located within the following region: PMSVSNDGKPSDLPHTFEYVVGHTMERSPMHVRNVGNPSDLPRTFEFMKG
Formulation:	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Concentration:	lot specific
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Predicted Protein Size:	69 kDa
Gene Name:	zinc finger protein 440
Database Link:	<a href="#">NP_689570</a> <a href="#">Entrez Gene 126070 Human</a> <a href="#">Q8IYI8</a>
Background:	May be involved in transcriptional regulation.
Synonyms:	FLJ37933
Note:	Immunogen Sequence Homology: Human: 100%
Protein Families:	Transcription Factors



[View online »](#)

**Product images:**


WB Suggested Anti-ZNF440 Antibody Titration:  
0.2-1 ug/ml; ELISA Titer: 1:312500; Positive  
Control: HepG2 cell lysate