

Product datasheet for **TA341502**


ZNF75A Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	WB
Recommended Dilution:	WB
Reactivity:	Human
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	The immunogen for anti-ZNF75A antibody: synthetic peptide directed towards the N terminal of human ZNF75A. Synthetic peptide located within the following region: FVLPKPKVISCLEQGEEPWVQVSPEFKDSAGKSPTGLKLNKNDTENHQPV
Formulation:	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose. <i>Note that this product is shipped as lyophilized powder to China customers.</i>
Concentration:	lot specific
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Predicted Protein Size:	35 kDa
Gene Name:	zinc finger protein 75a
Database Link:	NP_694573 Entrez Gene 7627 Human Q96N20
Background:	May be involved in transcriptional regulation.
Synonyms:	FLJ31529; MGC59959
Note:	Immunogen Sequence Homology: Human: 100%
Protein Families:	Transcription Factors



[View online »](#)

Product images:

WB Suggested Anti-ZNF75A Antibody Titration:
0.2-1 ug/ml; ELISA Titer: 1:12500; Positive Control:
MCF7 cell lysate