

Product datasheet for **TA341515**

ZNF100 Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	WB
Recommended Dilution:	WB
Reactivity:	Human
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	The immunogen for anti-ZNF100 antibody: synthetic peptide directed towards the N terminal of human ZNF100. Synthetic peptide located within the following region: MDDPRYGMCPKLGASGCPGAERSLLVQSYFEKGPLTFRDVAIEFSLEEWQ
Formulation:	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Concentration:	lot specific
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Predicted Protein Size:	63 kDa
Gene Name:	zinc finger protein 100
Database Link:	NP_775802 Entrez Gene 163227 Human Q8IYN0
Background:	May be involved in transcriptional regulation.
Synonyms:	FLJ44587
Note:	Immunogen Sequence Homology: Human: 100%
Protein Families:	Transcription Factors


[View online »](#)

Product images:


WB Suggested Anti-ZNF100 Antibody Titration:
0.2-1 ug/ml; Positive Control: Human Muscle