

## Product datasheet for **TA341516**


### ZNF596 Rabbit Polyclonal Antibody

#### Product data:

Product Type:	Primary Antibodies
Applications:	WB
Recommended Dilution:	WB
Reactivity:	Human
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	The immunogen for anti-ZNF596 antibody: synthetic peptide directed towards the C terminal of human ZNF596. Synthetic peptide located within the following region: AFNRSYNFRLHRRVHTGEKPYVCPLCGKAFSKFFNLRQHERHTTKAMNM
Formulation:	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose. <i>Note that this product is shipped as lyophilized powder to China customers.</i>
Concentration:	lot specific
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Predicted Protein Size:	58 kDa
Gene Name:	zinc finger protein 596
Database Link:	<a href="#">NP_775810</a> <a href="#">Entrez Gene 169270 Human</a> <a href="#">Q8TC21</a>
Background:	May be involved in transcriptional regulation.
Synonyms:	FLJ36123
Note:	Immunogen Sequence Homology: Dog: 100%; Pig: 100%; Goat: 100%; Horse: 100%; Human: 100%; Bovine: 100%; Zebrafish: 100%; Rabbit: 90%
Protein Families:	Druggable Genome, Transcription Factors



[View online »](#)

**Product images:**

WB Suggested Anti-ZNF596 Antibody Titration:  
0.2-1 ug/ml; ELISA Titer: 1:312500; Positive  
Control: Hela cell lysate