

Product datasheet for **TA341534**

ZNF735 Rabbit Polyclonal Antibody


Product data:

Product Type:	Primary Antibodies
Applications:	WB
Recommended Dilution:	WB
Reactivity:	Human
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	The immunogen for Anti-ZNF735 antibody is: synthetic peptide directed towards the middle region of Human ZNF735. Synthetic peptide located within the following region: NTQNKIFQTHKCVKVFSEKFSNSNRHNARYTGKKHLKCKKYGKSFCEMFSHL
Formulation:	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose. <i>Note that this product is shipped as lyophilized powder to China customers.</i>
Concentration:	lot specific
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Predicted Protein Size:	45 kDa
Gene Name:	zinc finger protein 735
Database Link:	NP_001152996 Entrez Gene 730291 Human P0CB33
Background:	May be involved in transcriptional regulation.
Synonyms:	ZNF735P
Note:	Immunogen Sequence Homology: Human: 100%



[View online »](#)

Product images:



Host: Rabbit
Target Name: ZNF735
Sample Tissue: Uterus Tumor Lysate
Antibody Dilution: 1.0 μ g/ml

Host: Rabbit; Target Name: ZNF735; Sample Tissue: Uterus Tumor lysates; Antibody Dilution: 1.0 ug/ml