

Product datasheet for **TA341541**


ZNF713 Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	WB
Recommended Dilution:	WB
Reactivity:	Human
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	The immunogen for anti-ZNF713 antibody: synthetic peptide directed towards the N terminal of human ZNF713. Synthetic peptide located within the following region: CKPEVIAQLELEEEWVIERDSLDDTHPDGENRPEIKKSTTSQNISDENQT
Formulation:	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Concentration:	lot specific
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Predicted Protein Size:	50 kDa
Gene Name:	zinc finger protein 713
Database Link:	NP_872439 Entrez Gene 349075 Human
Background:	May be involved in transcriptional regulation.
Synonyms:	FLJ39963
Note:	Immunogen Sequence Homology: Dog: 100%; Human: 100%; Bovine: 100%; Horse: 92%



[View online »](#)

Product images:

WB Suggested Anti-ZNF713 Antibody Titration:
0.2-1 ug/ml; ELISA Titer: 1:62500; Positive Control:
OVCAR-3 cell lysate