

Product datasheet for **TA341544**

ZNF17 Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	WB
Recommended Dilution:	WB
Reactivity:	Human
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	The immunogen for anti-ZNF17 antibody: synthetic peptide directed towards the N terminal of human ZNF17. Synthetic peptide located within the following region: NLTCMQGGKDFGSDSLQQALHSGWKPHRDTHGVEAFQSGQNNYSCTQC
Formulation:	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose. <i>Note that this product is shipped as lyophilized powder to China customers.</i>
Concentration:	lot specific
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Predicted Protein Size:	53 kDa
Gene Name:	zinc finger protein 17
Database Link:	NP_008890 Entrez Gene 7565 Human P17021
Background:	May be involved in transcriptional regulation.
Synonyms:	HPF3; KOX10
Note:	Immunogen Sequence Homology: Human: 100%



[View online »](#)

Product images:



WB Suggested Anti-ZNF17 Antibody Titration: 0.2-1 ug/ml; ELISA Titer: 1:12500; Positive Control: Human Small Intestine