

## Product datasheet for **TA341961**

### TMEM24 (C2CD2L) Rabbit Polyclonal Antibody

#### Product data:

Product Type:	Primary Antibodies
Applications:	WB
Recommended Dilution:	WB
Reactivity:	Human
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	The immunogen for anti-TMEM24 antibody: synthetic peptide directed towards the C terminal of human TMEM24. Synthetic peptide located within the following region: AGLSQSHDDLSNATATPSVRKKAGSFSRRLIKRFSFKSKPKANGNPSPQL
Formulation:	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose. <i>Note that this product is shipped as lyophilized powder to China customers.</i>
Concentration:	lot specific
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Predicted Protein Size:	76 kDa
Gene Name:	C2CD2 like
Database Link:	<a href="#">NP_055622</a> <a href="#">Entrez Gene 9854 Human</a> <a href="#">O14523</a>
Background:	The function of GADD45B is involved inThe regulation of growth and apoptosis.
Synonyms:	DLNB23; TMEM24
Note:	Immunogen Sequence Homology: Dog: 100%; Pig: 100%; Rat: 100%; Horse: 100%; Human: 100%; Mouse: 100%; Bovine: 100%; Rabbit: 100%; Guinea pig: 100%
Protein Families:	Transmembrane



[View online »](#)

## Product images:



WB Suggested Anti-TMEM24 Antibody Titration:  
0.2-1 ug/ml; Positive Control: Transfected 293T