

Product datasheet for **TA341973**

SUN2 Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	WB
Recommended Dilution:	WB
Reactivity:	Human
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	The immunogen for anti-UNC84B antibody: synthetic peptide directed towards the C terminal of human UNC84B. Synthetic peptide located within the following region: CWAFAQGPQGFAVVRLSARIRPTAVTLEHVPKALSPNSTISSAPKDFAIHG
Formulation:	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose. <i>Note that this product is shipped as lyophilized powder to China customers.</i>
Concentration:	lot specific
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Predicted Protein Size:	48 kDa
Gene Name:	Sad1 and UNC84 domain containing 2
Database Link:	NP_001186508 Entrez Gene 25777 Human Q9UH99



[View online »](#)

Background:

SUN1 (MIM 607723) and SUN2 are inner nuclear membrane (INM) proteins that play a major role in nuclear-cytoplasmic connection by formation of a 'bridge' across the nuclear envelope, known as the LINC complex, via interaction with the conserved luminal KASH domain of nesprins (e.g., SYNE1; MIM 608441) located in the outer nuclear membrane (ONM). The LINC complex provides a direct connection between the nuclear lamina and the cytoskeleton, which contributes to nuclear positioning and cellular rigidity (summary by Haque et al., 2010 [PubMed 19933576]). [supplied by OMIM, Nov 2010]. Transcript Variant: This variant (1) represents the longest transcript and encodes the longer isoform (a). Publication Note: This RefSeq record includes a subset of the publications that are available for this gene. Please see the Gene record to access additional publications. ##Evidence-Data-START## Transcript exon combination :: BC094797.1 [ECO:0000332] RNAseq introns :: single sample supports all introns ERS025087, ERS025093 [ECO:0000348] ##Evidence-Data-END## COMPLETENESS: complete on the 3' end.

Synonyms:


UNC84B

Note:

Immunogen Sequence Homology: Pig: 100%; Rat: 100%; Human: 100%; Mouse: 100%; Rabbit: 100%; Guinea pig: 100%; Dog: 93%; Horse: 93%; Bovine: 93%

Protein Families:

Transmembrane

Product images:

WB Suggested Anti-UNC84B Antibody Titration: 0.2-1 ug/ml; Positive Control: HepG2 cell lysate. SUN2 is supported by BioGPS gene expression data to be expressed in HepG2