

Product datasheet for **TA342106**

ZNF676 Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	WB
Recommended Dilution:	WB
Reactivity:	Human
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	The immunogen for anti-ZNF676 antibody: synthetic peptide directed towards the middle region of human ZNF676. Synthetic peptide located within the following region: KPYKCEECGKGFSSVSTLNTHKAIHAEKPYKCEECGKASNSSSKLMEHK
Formulation:	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Concentration:	lot specific
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Predicted Protein Size:	67 kDa
Gene Name:	zinc finger protein 676
Database Link:	NP_001001411 Entrez Gene 163223 Human Q8N7Q3
Background:	ZNF676 may be involved in transcriptional regulation.
Synonyms:	zinc finger protein 676



Note: Immunogen Sequence Homology: Human: 100%

Product images:



WB Suggested Anti-ZNF676 Antibody Titration:
0.2-1 ug/ml; ELISA Titer: 1: 1562500; Positive
Control: Jurkat cell lysate