

Product datasheet for **TA342108**

ZNF782 Rabbit Polyclonal Antibody


Product data:

Product Type:	Primary Antibodies
Applications:	WB
Recommended Dilution:	WB
Reactivity:	Human
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	The immunogen for anti-ZNF782 antibody: synthetic peptide directed towards the N terminal of human ZNF782. Synthetic peptide located within the following region: TNKLLTTEQEISGKPHNRDINIFRARMMPCKCDIAGSACQGLSLMAPHCQ
Formulation:	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose. <i>Note that this product is shipped as lyophilized powder to China customers.</i>
Concentration:	lot specific
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Predicted Protein Size:	81 kDa
Gene Name:	zinc finger protein 782
Database Link:	NP_001001662 Entrez Gene 158431 Human Q6ZMW2
Background:	ZNF782 may be involved in transcriptional regulation.
Synonyms:	FLJ16636
Note:	Immunogen Sequence Homology: Human: 100%; Dog: 79%
Protein Families:	Transcription Factors



[View online »](#)

Product images:



WB Suggested Anti-ZNF782 Antibody Titration:
0.2-1 ug/ml; ELISA Titer: 1: 62500; Positive
Control: Jurkat cell lysate