


## Product datasheet for **TA342121**

### ZNF599 Rabbit Polyclonal Antibody

#### Product data:

Product Type:	Primary Antibodies
Applications:	WB
Recommended Dilution:	WB
Reactivity:	Human
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	The immunogen for anti-ZNF599 antibody: synthetic peptide directed towards the N terminal of human ZNF599. Synthetic peptide located within the following region: MAAPALALVSFEDVVVFTTGEEWGHLDLAQRTLYQEVMLETCRLVSLGH
Formulation:	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Concentration:	lot specific
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Predicted Protein Size:	67 kDa
Gene Name:	zinc finger protein 599
Database Link:	<a href="#">NP_001007249</a> <a href="#">Entrez Gene 148103 Human</a> <a href="#">Q96NL3</a>
Background:	ZNF599 may be involved in transcriptional regulation.
Synonyms:	FLJ50920
Note:	Immunogen Sequence Homology: Human: 100%; Bovine: 92%


[View online »](#)

**Product images:**

WB Suggested Anti-ZNF599 Antibody Titration:  
0.2-1 ug/ml; ELISA Titer: 1: 62500; Positive  
Control: 293T cell lysate