

Product datasheet for **TA342126**

ZNF766 Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	WB
Recommended Dilution:	WB
Reactivity:	Human
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	The immunogen for anti-ZNF766 antibody: synthetic peptide directed towards the N terminal of human ZNF766. Synthetic peptide located within the following region: MKVAKNPDRWEGIKDINTGRSCAVRSKAGNKPITNQLGLTFQLPLPELEI
Formulation:	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Concentration:	lot specific
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Predicted Protein Size:	54 kDa
Gene Name:	zinc finger protein 766
Database Link:	NP_001010851 Entrez Gene 90321 Human Q5HY98
Background:	ZNF766 belongs to the krueppel C2H2-type zinc-finger protein family. It contains 10 C2H2-type zinc fingers and 1 KRAB domain. ZNF766 may be involved in transcriptional regulation.
Synonyms:	zinc finger protein 766
Note:	Immunogen Sequence Homology: Human: 100%


[View online »](#)

Product images:


WB Suggested Anti-ZNF766 Antibody Titration:
0.2-1 ug/ml; Positive Control: Hela cell lysate