

Product datasheet for **TA342318**

Zscan12 Rabbit Polyclonal Antibody


Product data:

Product Type:	Primary Antibodies
Applications:	WB
Recommended Dilution:	WB
Reactivity:	Mouse
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	The immunogen for anti-Zscan12 antibody: synthetic peptide corresponding to a region of Mouse. Synthetic peptide located within the following region: NSSLATHQETHHKEKPFTQSGPIQQQRNHTKEKPYKCSVCGKAFIQKISL
Formulation:	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose. <i>Note that this product is shipped as lyophilized powder to China customers.</i>
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Predicted Protein Size:	57 kDa
Gene Name:	zinc finger and SCAN domain containing 12
Database Link:	NP_057893 Entrez Gene 22758 Mouse Q9Z1D7
Background:	May be involved in transcriptional regulation.
Synonyms:	dj29K1.2; KIAA0426; OTTHUMP00000016205; ZFP96; ZNF29K1; ZNF96; ZNF305
Note:	Immunogen Sequence Homology: Mouse: 100%; Rat: 85%



[View online »](#)

Product images:



WB Suggested Anti-Zscan12 Antibody Titration:
0.2-1 ug/ml; Positive Control: Mouse Thymus