

Product datasheet for **TA342409**

ZNF96 (ZSCAN12) Rabbit Polyclonal Antibody


Product data:

Product Type:	Primary Antibodies
Applications:	WB
Recommended Dilution:	WB
Reactivity:	Human
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	The immunogen for anti-ZSCAN12 antibody: synthetic peptide directed towards the N terminal of human ZSCAN12. Synthetic peptide located within the following region: QEMFLQETVRLRKEGEPMSLSLQSMKAQPKYESPELESQQEQVLDVETGNE
Formulation:	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose. <i>Note that this product is shipped as lyophilized powder to China customers.</i>
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Predicted Protein Size:	66 kDa
Gene Name:	zinc finger and SCAN domain containing 12
Database Link:	NP_001034732 Entrez Gene 9753 Human
Background:	ZSCAN12 may be involved in transcriptional regulation.
Synonyms:	dj29K1.2; KIAA0426; ZFP96; zinc finger and SCAN domain containing 12; zinc finger protein 96; zinc finger protein 305; ZNF29K1; ZNF96; ZNF305; ZSCAN12
Note:	Immunogen Sequence Homology: Human: 100%; Pig: 93%; Bovine: 77%
Protein Families:	Transcription Factors



[View online »](#)

Product images:



WB Suggested Anti-ZSCAN12 Antibody Titration:
5.0 ug/ml; Positive Control: Jurkat cell lysate