

Product datasheet for **TA342423**

ZNF91 Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	WB
Recommended Dilution:	WB
Reactivity:	Human
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	The immunogen for anti-ZNF91 antibody: synthetic peptide directed towards the N terminal of human ZNF91. Synthetic peptide located within the following region: HKIIHAGEKLYKCEECGKAFNRSSNLTIHKFIHTGEKPYKCEECGKAFNW
Formulation:	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose. <i>Note that this product is shipped as lyophilized powder to China customers.</i>
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Predicted Protein Size:	137 kDa
Gene Name:	zinc finger protein 91
Database Link:	NP_003421 Entrez Gene 7644 Human Q05481



[View online »](#)

Background:

The ZNF91 gene encodes a zinc finger protein of the KRAB (Kruppel-associated box) subfamily (Bellefroid et al., 1991, 1993 [PubMed 2023909] [PubMed 8467795]). [supplied by OMIM, May 2010]. Transcript Variant: This variant (1) represents the longer transcript and encodes the longer isoform (1). Sequence Note: This RefSeq record was created from transcript and genomic sequence data to make the sequence consistent with the reference genome assembly. The genomic coordinates used for the transcript record were based on transcript alignments. ##Evidence-Data-START## Transcript exon combination :: DB068854.1, AK291256.1 [ECO:0000332] ##Evidence-Data-END## COMPLETENESS: complete on the 3' end.

Synonyms:


HPF7; HTF10

Note:

Immunogen Sequence Homology: Human: 100%; Rat: 93%; Mouse: 93%

Protein Families:

Transcription Factors

Product images:

WB Suggested Anti-ZNF91 Antibody Titration: 0.2-1 ug/ml; ELISA Titer: 1: 312500; Positive Control: Human Muscle