

Product datasheet for **TA342425**

ZNF131 Rabbit Polyclonal Antibody


Product data:

Product Type:	Primary Antibodies
Applications:	WB
Recommended Dilution:	WB
Reactivity:	Human
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	The immunogen for anti-ZNF131 antibody: synthetic peptide directed towards the middle region of human ZNF131. Synthetic peptide located within the following region: DLLQDSQVHDSHMSSELPEQVQVSYLEVGRIQTEEGTEVHVEELHVERVNQ
Formulation:	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose. <i>Note that this product is shipped as lyophilized powder to China customers.</i>
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Predicted Protein Size:	67 kDa
Gene Name:	zinc finger protein 131
Database Link:	NP_003423 Entrez Gene 7690 Human P52739
Background:	ZNF131 may be involved in transcriptional regulation. It plays a role during development and organogenesis as well as in the function of the adult central nervous system.
Synonyms:	pHZ-10; ZBTB35
Note:	Immunogen Sequence Homology: Dog: 100%; Pig: 100%; Rat: 100%; Horse: 100%; Human: 100%; Mouse: 100%; Bovine: 100%; Rabbit: 100%; Guinea pig: 100%
Protein Families:	Transcription Factors



[View online »](#)

Product images:



WB Suggested Anti-ZNF131 Antibody Titration:
0.2-1 ug/ml; ELISA Titer: 1: 62500; Positive
Control: Human Small Intestine