

Product datasheet for **TA342429**

ZKSCAN1 Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	WB
Recommended Dilution:	WB
Reactivity:	Human
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	The immunogen for anti-ZKSCAN1 antibody: synthetic peptide directed towards the N terminal of human ZKSCAN1. Synthetic peptide located within the following region: NTKEQILELLVLEQFLSILPKELQVWLQEYRPDSGEEAVTLLEDLELDLS
Formulation:	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose. <i>Note that this product is shipped as lyophilized powder to China customers.</i>
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Predicted Protein Size:	63 kDa
Gene Name:	zinc finger with KRAB and SCAN domains 1
Database Link:	NP_003430 Entrez Gene 7586 Human P17029



[View online »](#)

Background:

The ZKSCAN1 gene encodes a transcriptional regulator of the KRAB (Kruppel-associated box) subfamily of zinc finger proteins, which contain repeated Cys2-His2 (C2H2) zinc finger domains that are connected by conserved sequences, called H/C links (summarized by Tommerup and Vissing, 1995 [PubMed 7557990]). Transcriptional regulatory proteins containing tandemly repeated zinc finger domains are thought to be involved in both normal and abnormal cellular proliferation and differentiation. See ZNF91 (MIM 603971) for general information on zinc finger proteins. [supplied by OMIM, Jul 2010]. Transcript Variant: This variant (3) differs in its 5' UTR, lacks a portion of the 5' coding region and initiates translation at a downstream in-frame start codon, compared to variant 1. The encoded isoform (c) is shorter at the N-terminus, compared to isoform a. Sequence Note: This RefSeq record was created from transcript and genomic sequence data to make the sequence consistent with the reference genome assembly. The genomic coordinates used for the transcript record were based on transcript alignments. ##Evidence-Data-START## Transcript exon combination :: AK091471.1, DA499270.1 [ECO:0000332] RNAseq introns :: single sample supports all introns ERS025081 [ECO:0000348] ##Evidence-Data-END## COMPLETENESS: complete on the 3' end.

Synonyms:


9130423L19Rik; KOX18; PHZ-37; ZNF36; ZNF139; ZSCAN33

Note:

Immunogen Sequence Homology: Dog: 100%; Pig: 100%; Rat: 100%; Horse: 100%; Human: 100%; Mouse: 100%; Bovine: 100%; Rabbit: 100%; Guinea pig: 100%

Protein Families:

Transcription Factors

Product images:

WB Suggested Anti-ZKSCAN1 Antibody Titration: 5.0 ug/ml; Positive Control: Jurkat cell lysate