

Product datasheet for TA342499

GLYCTK Rabbit Polyclonal Antibody

Product data:

Product Type: Primary Antibodies

Applications: WB

Recommended Dilution: WB

Reactivity: Human

Host: Rabbit

Isotype: IgG

Clonality: Polyclonal

Immunogen: The immunogen for anti-GLYCTK antibody: synthetic peptide directed towards the N terminal

of human GLYCTK. Synthetic peptide located within the following region:

IQQLAEGLTADDLLLVLISGGGSALLPAPIPPVTLEEKQTLTRLLAARGA

Formulation: Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2%

sucrose.

Note that this product is shipped as lyophilized powder to China customers.

Conjugation: Unconjugated

Storage: Store at -20°C as received.

Stability: Stable for 12 months from date of receipt.

Predicted Protein Size: 55 kDa

Gene Name: glycerate kinase

Database Link: NP 660305

Entrez Gene 132158 Human

Q8IVS8

Background: This locus encodes a member of the glycerate kinase type-2 family. The encoded enzyme

catalyzes the phosphorylation of (R)-glycerate and may be involved in serine degradation and fructose metabolism. Decreased activity of the encoded enzyme may be associated with the disease D-glyceric aciduria. Alternatively spliced transcript variants have been described.

[provided by RefSeq, Jan 2009]

Synonyms: HBeAgBP4A; HBEBP2; HBEBP4



OriGene Technologies, Inc. 9620 Medical Center Drive, Ste 200

CN: techsupport@origene.cn

Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com



GLYCTK Rabbit Polyclonal Antibody - TA342499

Note: Immunogen Sequence Homology: Dog: 100%; Pig: 100%; Rat: 100%; Horse: 100%; Human:


100%; Mouse: 100%; Bovine: 100%; Rabbit: 100%; Guinea pig: 100%

Protein Families: Transcription Factors

Protein Pathways: Glycerolipid metabolism, Glycine, serine and threonine metabolism, Glyoxylate and

dicarboxylate metabolism, Metabolic pathways

Product images:



WB Suggested Anti-GLYCTK Antibody Titration: 1.25 ug/ml; Positive Control: HepG2 cell lysate