

Product datasheet for **TA343194**

PITPNB Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	WB
Recommended Dilution:	WB
Reactivity:	Human
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	The immunogen for anti-PITPNB antibody is: synthetic peptide directed towards the N-terminal region of Human PITPNB. Synthetic peptide located within the following region: EGIEVLKNEPYEKDGEKQYTHKIYHLKSKVPAFVRMIAPEGSLVFHEKA
Formulation:	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose. <i>Note that this product is shipped as lyophilized powder to China customers.</i>
Purification:	Affinity Purified
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Predicted Protein Size:	30 kDa
Gene Name:	phosphatidylinositol transfer protein beta
Database Link:	NP_036531 Entrez Gene 23760 Human P48739
Background:	The function of this protein remains unknown.
Synonyms:	PI-TP-beta; PtdInsTP; VIB1B
Note:	Immunogen Sequence Homology: Dog: 100%; Pig: 100%; Rat: 100%; Horse: 100%; Human: 100%; Bovine: 100%; Rabbit: 100%; Zebrafish: 100%; Mouse: 93%; Guinea pig: 93%



[View online »](#)

Product images:

Host: Rabbit
Target Name: PITPNB
Sample Tissue: Fetal Brain Lysate
Antibody Dilution: 1.0 μ g/ml

Host: Rabbit; Target Name: PITPNB; Sample Tissue: Fetal Brain lysates; Antibody Dilution: 1.0 ug/ml