

Product datasheet for **TA343405**


ZNF189 Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	WB
Recommended Dilution:	WB
Reactivity:	Human
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	The immunogen for anti-ZNF189 antibody: synthetic peptide directed towards the N terminal of human ZNF189. Synthetic peptide located within the following region: GLLTFEDVAVFFTQEWDYLDPAQRSLYKDVM MENYGNLVS LDVLRDKD
Formulation:	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose. <i>Note that this product is shipped as lyophilized powder to China customers.</i>
Purification:	Affinity Purified
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Predicted Protein Size:	73 kDa
Gene Name:	zinc finger protein 189
Database Link:	NP_003443 Entrez Gene 7743 Human O75820
Background:	ZNF189 may be involved in transcriptional regulation.
Synonyms:	OTTHUMP00000021808; OTTHUMP00000021810; zinc finger protein 189
Note:	Immunogen Sequence Homology: Dog: 100%; Rat: 100%; Horse: 100%; Human: 100%; Bovine: 100%; Pig: 92%
Protein Families:	Transcription Factors



[View online »](#)

Product images:

WB Suggested Anti-ZNF189 Antibody Titration:
0.2-1 ug/ml; ELISA Titer: 1: 1562500; Positive
Control: HepG2 cell lysate; ZNF189 is supported
by BioGPS gene expression data to be expressed
in HepG2