

Product datasheet for **TA343407**

Zfp202 Rabbit Polyclonal Antibody


Product data:

Product Type:	Primary Antibodies
Applications:	WB
Recommended Dilution:	WB
Reactivity:	Mouse
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	The immunogen for Anti-Zfp202 antibody is: synthetic peptide directed towards the N-terminal region of Mouse Zfp202. Synthetic peptide located within the following region: MATALEPEDQDLWEEEGVLMVKLEDDFTCGPESALEGDDPVLETSHQNFR
Formulation:	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Purification:	Affinity Purified
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Predicted Protein Size:	30 kDa
Gene Name:	zinc finger protein 202
Database Link:	NP_109638 Entrez Gene 80902 Mouse
Background:	The function of this protein remains unknown.
Synonyms:	C130037E22Rik; RGD1562294; Znf202



Note: Immunogen Sequence Homology: Human: 100%; Dog: 93%; Pig: 93%; Horse: 93%; Mouse: 93%; Bovine: 93%; Guinea pig: 93%; Rat: 89%; Rabbit: 86%

Product images:



Host: Rabbit; Target Name: Zfp202; Sample
Tissue: Mouse Kidney lysates; Antibody
Dilution: 1.0 ug/ml