

Product datasheet for **TA343414**

Klf11 Rabbit Polyclonal Antibody


Product data:

Product Type:	Primary Antibodies
Applications:	WB
Recommended Dilution:	WB
Reactivity:	Rat
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	The immunogen for Anti-Klf11 antibody is: synthetic peptide directed towards the N-terminal region of Rat Klf11. Synthetic peptide located within the following region: AADIMDICESILERKRHDSERSTCSILEQTDIEAVEALVCMSSWGQRSQV
Formulation:	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose. <i>Note that this product is shipped as lyophilized powder to China customers.</i>
Purification:	Affinity Purified
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Predicted Protein Size:	55 kDa
Gene Name:	Kruppel-like factor 11
Database Link:	NP_001032431 Entrez Gene 313994 Rat
Background:	The function of this protein remains unknown.
Synonyms:	FKLF; FKLF1; MODY7; TIEG-2; TIEG2; Tieg3
Note:	Immunogen Sequence Homology: Dog: 100%; Pig: 100%; Rat: 100%; Horse: 100%; Human: 100%; Mouse: 100%; Bovine: 100%; Rabbit: 100%; Guinea pig: 92%



[View online »](#)

Product images:



WB Suggested Anti-Klf11 Antibody; Titration: 1.0 ug/ml; Positive Control: Rat Stomach