


Product datasheet for **TA343601**

ZNF193 (ZSCAN9) Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	WB
Recommended Dilution:	WB
Reactivity:	Human
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	The immunogen for anti-ZNF193 antibody: synthetic peptide directed towards the N terminal of human ZNF193. Synthetic peptide located within the following region: MNTNSKEVLSLGVQVPEAWELLTMKVEAKSHLQWQESRLKRSNPLAREI
Formulation:	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Purification:	Affinity Purified
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Predicted Protein Size:	46 kDa
Gene Name:	zinc finger and SCAN domain containing 9
Database Link:	NP_006290 Entrez Gene 7746 Human O15535
Background:	ZNF193 belongs to the krueppel C2H2-type zinc-finger protein family. It contains 5 C2H2-type zinc fingers and 1 SCAN box domain. ZNF193 may be involved in transcriptional regulation.
Synonyms:	PRD51; ZNF193
Note:	Immunogen Sequence Homology: Human: 100%
Protein Families:	Transcription Factors


[View online »](#)

Product images:


WB Suggested Anti-ZNF193 Antibody Titration:
 0.2-1 ug/ml; ELISA Titer: 1: 312500; Positive
 Control: HT1080 cell lysate