

Product datasheet for **TA343630**

iASPP (PPP1R13L) Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	WB
Recommended Dilution:	WB
Reactivity:	Human
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	The immunogen for anti-PPP1R13L antibody: synthetic peptide directed towards the middle region of human PPP1R13L. Synthetic peptide located within the following region: GLMNSGAVYALWDYSAEFGDELSFREGESVTVLRRDGPETDWWWAALHG
Formulation:	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose. <i>Note that this product is shipped as lyophilized powder to China customers.</i>
Purification:	Affinity Purified
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Predicted Protein Size:	89 kDa
Gene Name:	protein phosphatase 1 regulatory subunit 13 like
Database Link:	NP_006654 Entrez Gene 10848 Human Q8WUF5
Background:	IASPP is one of the most evolutionarily conserved inhibitors of p53 (TP53), whereas ASPP1 and ASPP2 are activators of p53. IASPP is one of the most evolutionarily conserved inhibitors of p53 (TP53; MIM 191170), whereas ASPP1 (MIM 606455) and ASPP2 (MIM 602143) are activators of p53. [supplied by OMIM]. Publication Note: This RefSeq record includes a subset of the publications that are available for this gene. Please see the Entrez Gene record to access additional publications.



[View online »](#)

Synonyms: IASPP; NKIP1; RAI; RAI4

Note: Immunogen Sequence Homology: Dog: 100%; Rat: 100%; Horse: 100%; Human: 100%; Mouse: 100%; Rabbit: 100%; Pig: 93%; Bovine: 93%; Guinea pig: 93%

Protein Families: Transcription Factors

Product images:



WB Suggested Anti-PPP1R13L Antibody Titration:
0.2-1 ug/ml; ELISA Titer: 1: 312500; Positive
Control: Human Lung