

Product datasheet for **TA343874**

KLHDC8A Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	WB
Recommended Dilution:	WB
Reactivity:	Human
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	The immunogen for anti-KLHDC8A antibody: synthetic peptide directed towards the C terminal of human KLHDC8A. Synthetic peptide located within the following region: PGKKNKWEILPAMPTPRCACSSIVVKNCLLAVGGVNQGLSDAVEALCVSDS
Formulation:	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Purification:	Affinity Purified
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Predicted Protein Size:	39 kDa
Gene Name:	kelch domain containing 8A
Database Link:	NP_060673 Entrez Gene 55220 Human Q8IYD2
Background:	The function remains unknown.
Synonyms:	FLJ10748
Note:	Immunogen Sequence Homology: Pig: 100%; Rat: 100%; Human: 100%; Bovine: 100%; Rabbit: 100%; Guinea pig: 100%; Dog: 93%; Horse: 93%; Mouse: 93%



[View online »](#)

Product images:



WB Suggested Anti-KLHDC8A Antibody Titration:
1.25 ug/ml; Positive Control: Jurkat cell lysate