

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

Product datasheet for TA343903

SF4 (SUGP1) Rabbit Polyclonal Antibody


Product data:

Product Type:	Primary Antibodies
Applications:	WB
Recommended Dilution:	WB
Reactivity:	Human
Host:	Rabbit
lsotype:	IgG
Clonality:	Polyclonal
Immunogen:	The immunogen for anti-SF4 antibody: synthetic peptide directed towards the N terminal of human SF4. Synthetic peptide located within the following region: MSLKMDNRDVAGKANRWFGVAPPKSGKMNMNILHQEELIAQKKREIEAKM
Formulation:	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose. Note that this product is shipped as lyophilized powder to China customers.
Purification:	Affinity Purified
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Predicted Protein Size:	72 kDa
Gene Name:	SURP and G-patch domain containing 1
Database Link:	<u>NP_757386</u> <u>Entrez Gene 57794 Human</u> <u>Q8IWZ8</u>
Background:	SF4 is a member of the SURP family of splicing factors.SF4 is a member of the SURP family of splicing factors. [supplied by OMIM]
Synonyms:	F23858; RBP; SF4
Note:	lmmunogen Sequence Homology: Dog: 100%; Pig: 100%; Rat: 100%; Horse: 100%; Human: 100%; Mouse: 100%; Bovine: 100%; Zebrafish: 100%; Guinea pig: 100%; Yeast: 83%



This product is to be used for laboratory only. Not for diagnostic or therapeutic use. ©2022 OriGene Technologies, Inc., 9620 Medical Center Drive, Ste 200, Rockville, MD 20850, US

Product images:



WB Suggested Anti-SF4 Antibody Titration: 0.2-1 ug/ml; Positive Control: HepG2 cell lysate

This product is to be used for laboratory only. Not for diagnostic or therapeutic use. ©2022 OriGene Technologies, Inc., 9620 Medical Center Drive, Ste 200, Rockville, MD 20850, US