

## Product datasheet for **TA343910**

### NXF3 Rabbit Polyclonal Antibody

#### Product data:

Product Type:	Primary Antibodies
Applications:	WB
Recommended Dilution:	WB
Reactivity:	Human
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	The immunogen for anti-NXF3 antibody: synthetic peptide directed towards the N terminal of human NXF3. Synthetic peptide located within the following region: SLPSGHTTGHTDQVVQRRARCWDIYQRRFSSRSEPVNPGMHSSSHQQQDG
Formulation:	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Purification:	Affinity Purified
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Predicted Protein Size:	60 kDa
Gene Name:	nuclear RNA export factor 3
Database Link:	<a href="#">NP_071335</a> <a href="#">Entrez Gene 56000 Human</a> <a href="#">Q9H4D5</a>


[View online »](#)

**Background:**


NXF3 is one member of a family of nuclear RNA export factors. Common domain features of this family are a noncanonical RNP-type RNA-binding domain (RBD), 4 leucine-rich repeats (LRRs), a nuclear transport factor 2 (NTF2)-like domain that allows heterodimerization with NTF2-related export protein-1 (NXT1), and a ubiquitin-associated domain that mediates interactions with nucleoporins. The LRRs and NTF2-like domains are required for export activity. NXF3 has shortened LRR and ubiquitin-associated domains and its RDB is unable to bind RNA. It is located in the nucleoplasm but is not associated with either the nuclear envelope or the nucleolus. This gene is one member of a family of nuclear RNA export factor genes. Common domain features of this family are a noncanonical RNP-type RNA-binding domain (RBD), 4 leucine-rich repeats (LRRs), a nuclear transport factor 2 (NTF2)-like domain that allows heterodimerization with NTF2-related export protein-1 (NXT1), and a ubiquitin-associated domain that mediates interactions with nucleoporins. The LRRs and NTF2-like domains are required for export activity. Alternative splicing seems to be a common mechanism in this gene family. The encoded protein of this gene has shortened LRR and ubiquitin-associated domains and its RDB is unable to bind RNA. It is located in the nucleoplasm but is not associated with either the nuclear envelope or the nucleolus.

**Synonyms:**


TAPL-3; TAPL3

**Note:**

Immunogen Sequence Homology: Human: 100%

**Product images:**


WB Suggested Anti-NXF3 Antibody Titration: 0.2-1 ug/ml; Positive Control: HepG2 cell lysate



Lanes: Lane 1: 7ug HeLa cytosolic lysate; Lane 2: 7ug HeLa nuclear lysate; Primary Antibody Dilution: 1: 300; Secondary Antibody: Anti-rabbit-HRP; Secondary Antibody Dilution: 1: 500; Gene Name: NXF3