

Product datasheet for TA344233

IMPA1 Rabbit Polyclonal Antibody

Product data:

Product Type: Primary Antibodies

Applications: WB

Recommended Dilution: WB

Reactivity: Human

Host: Rabbit

Isotype: IgG

Clonality: Polyclonal

Immunogen: The immunogen for anti-IMPA1 antibody: synthetic peptide directed towards the middle

region of human IMPA1. Synthetic peptide located within the following region:

IVTEAGGVLMDVTGGPFDLMSRRVIAANNRILAERIAKEIQVIPLQRDDE

Formulation: Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2%

sucrose.

Note that this product is shipped as lyophilized powder to China customers.

Purification: Affinity Purified

Conjugation: Unconjugated

Store at -20°C as received.

Stability: Stable for 12 months from date of receipt.

Predicted Protein Size: 30 kDa

Gene Name: inositol monophosphatase 1

Database Link: NP 005527

Entrez Gene 3612 Human

P29218

Background: IMPA1 is responsible for the provision of inositol required for synthesis of

phosphatidylinositol and polyphosphoinositides and has been implicated as the pharmacological target for lithium action in brain. IMPA1 can use myo-inositol

monophosphates, myo-inositol-1,3-diphosphate, myo-inositol-1,4-diphosphate, scyllo-inositol-

phosphate, glucose-1-phosphate, glucose-6-phosphate, fructose-1-phosphate, beta-

glycerophosphate, and 2'-AMP as substrates.



OriGene Technologies, Inc. 9620 Medical Center Drive, Ste 200

CN: techsupport@origene.cn

Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com





Synonyms: IMP; IMPA


Note: Immunogen Sequence Homology: Human: 100%; Dog: 93%; Pig: 93%; Bovine: 93%; Rat: 92%;

Guinea pig: 86%; Mouse: 85%; Horse: 79%; Rabbit: 79%

Protein Families: Druggable Genome

Protein Pathways: Inositol phosphate metabolism, Metabolic pathways, Phosphatidylinositol signaling system

Product images:



WB Suggested Anti-IMPA1 Antibody Titration: 0.2-1 ug/ml; ELISA Titer: 1:1562500; Positive Control: Human heart