

Product datasheet for TA344419

RBFOX1 Rabbit Polyclonal Antibody

Product data:

Product Type: Primary Antibodies

Applications: WB

Recommended Dilution: WB

Reactivity: Human

Host: Rabbit

Isotype: IgG

Clonality: Polyclonal

Immunogen: The immunogen for Anti-A2BP1 antibody is: synthetic peptide directed towards the N-

terminal region of Human A2BP1. Synthetic peptide located within the following region:

NCEREQLRGNQEAAAAPDTMAQPYASAQFAPPQNGIPAEYTAPHPHPAPE

Formulation: Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2%

sucrose.

Note that this product is shipped as lyophilized powder to China customers.

Purification: Affinity Purified

Conjugation: Unconjugated

Store at -20°C as received.

Stability: Stable for 12 months from date of receipt.

Predicted Protein Size: 40 kDa

Gene Name: RNA binding protein, fox-1 homolog 1

Database Link: NP 001135805

Entrez Gene 54715 Human

Q9NWB1



OriGene Technologies, Inc. 9620 Medical Center Drive, Ste 200

CN: techsupport@origene.cn

Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com




Background: Ataxin-2 binding protein 1(A2BP1) has an RNP motif that is highly conserved among RNA-

binding proteins. This protein binds to the C-terminus of ataxin-2 and may contribute to the restricted pathology of spinocerebellar ataxia type 2 (SCA2). Ataxin-2 is the gene product of the SCA2 gene which causes familial neurodegenerative diseases. Ataxin-2 binding protein 1 and ataxin-2 are both localized in the trans-Golgi network. Several alternatively spliced transcript variants encoding different isoforms have been found for this gene. Additional transcript variants have been found but their full length nature has not been determined.

Synonyms: 2BP1; A2BP1; FOX-1; FOX1; HRNBP1


Note: Immunogen Sequence Homology: Human: 100%; Mouse: 100%; Bovine: 100%

Product images:



WB Suggested Anti-A2BP1 Antibody Titration: 0.2-1 ug/ml; ELISA Titer: 1:62500; Positive Control: MCF7 cell lysate

A2BP1



2Human 293T; Host: Rabbit; Target Name: HIRIP3; Sample Tissue: 293T; Antibody Dilution: 1.0 ug/ml