

Product datasheet for TA344534

ARIH2 Rabbit Polyclonal Antibody

Product data:

Product Type: Primary Antibodies

Applications: WB

Recommended Dilution: WB

Reactivity: Rat, Human

Host: Rabbit

Isotype: IgG

Clonality: Polyclonal

Immunogen: The immunogen for Anti-Arih2 antibody is: synthetic peptide directed towards the C-terminal

region of Rat Arih2. Synthetic peptide located within the following region: YYMESGPRKKLFEYQQAQLEAEIENLSWKVERADSYDRGDLENQMHIAEQ

Formulation: Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2%

sucrose.

Note that this product is shipped as lyophilized powder to China customers.

Purification: Affinity Purified

Conjugation: Unconjugated

Store at -20°C as received.

Stability: Stable for 12 months from date of receipt.

Predicted Protein Size: 54 kDa

Gene Name: ariadne RBR E3 ubiquitin protein ligase 2

Database Link: NP 006312

Entrez Gene 316005 RatEntrez Gene 10425 Human

O95376

Background: The function of this protein remains unknown.

Synonyms: ARI2; TRIAD1

Note: Immunogen Sequence Homology: Dog: 100%; Pig: 100%; Rat: 100%; Horse: 100%; Human:

100%; Mouse: 100%; Bovine: 100%; Rabbit: 100%; Guinea pig: 100%

Protein Families: Druggable Genome




OriGene Technologies, Inc. 9620 Medical Center Drive, Ste 200

CN: techsupport@origene.cn

Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com



Product images:



Host: Rabbit; Target Name: Arih2; Sample Tissue: Rat Testis lysates; Antibody Dilution: 1.0 ug/ml