

Product datasheet for **TA344539**

THRAP5 (MED16) Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	WB
Recommended Dilution:	WB
Reactivity:	Human
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	The immunogen for anti-MED16 antibody: synthetic peptide directed towards the C terminal of human MED16. Synthetic peptide located within the following region: MSLLFRLTLKLWICCRDEGPASEPDEALVDECCLLPSQLLIPSLDWLPAS
Formulation:	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose. <i>Note that this product is shipped as lyophilized powder to China customers.</i>
Purification:	Affinity Purified
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Predicted Protein Size:	97 kDa
Gene Name:	mediator complex subunit 16
Database Link:	NP_005472 Entrez Gene 10025 Human Q9Y2X0
Background:	THRAP5 is part of the human thyroid hormone receptor-associated protein (TRAP)-Mediator family which acts as a coactivator for a broad range of nuclear hormone receptors as well as other classes of transcriptional activators.
Synonyms:	DRIP92; THRAP5; TRAP95
Note:	Immunogen Sequence Homology: Pig: 100%; Human: 100%; Guinea pig: 100%; Dog: 92%; Rat: 92%; Bovine: 92%; Mouse: 85%



[View online »](#)

Protein Families: Druggable Genome, Transcription Factors

Product images:



WB Suggested Anti-MED16 Antibody Titration:
0.2-1 ug/ml; ELISA Titer: 1:62500; Positive Control:
Human Placenta