

Product datasheet for **TA344573**

SPDYE3 Rabbit Polyclonal Antibody


Product data:

Product Type:	Primary Antibodies
Applications:	WB
Recommended Dilution:	WB
Reactivity:	Human
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	The immunogen for Anti-SPDYE3 antibody is: synthetic peptide directed towards the C-terminal region of Human SPDYE3. Synthetic peptide located within the following region: AETLCGLKMKAKRRRVSLVLPYYEAFNRLLEDPIKRFLLAWDKDLRVSD
Formulation:	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Purification:	Affinity Purified
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Predicted Protein Size:	18 kDa
Gene Name:	speedy/RINGO cell cycle regulator family member E3
Database Link:	NP_001004351 Entrez Gene 441272 Human A6NKU9
Background:	The function of this protein remains unknown.
Synonyms:	MGC57359; MGC119294
Note:	Immunogen Sequence Homology: Human: 100%; Dog: 91%; Pig: 91%; Rat: 91%; Mouse: 91%; Bovine: 91%; Rabbit: 91%; Guinea pig: 91%



[View online »](#)

Product images:



Host: Rabbit; Target Name: SPDYE3; Sample
Tissue: OVCAR-3 whole cell lysates; Antibody
Dilution: 1.0 ug/ml