

Product datasheet for **TA344792**

PCNP Rabbit Polyclonal Antibody


Product data:

Product Type:	Primary Antibodies
Applications:	WB
Recommended Dilution:	WB
Reactivity:	Human
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	The immunogen for anti-PCNP antibody: synthetic peptide directed towards the middle region of human PCNP. Synthetic peptide located within the following region: AKMRMKNIGRDTPTSAGPNSFNKGKHGFS DNQKLWERNIKSHLGNVHDQD
Formulation:	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose. <i>Note that this product is shipped as lyophilized powder to China customers.</i>
Purification:	Affinity Purified
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Predicted Protein Size:	19 kDa
Gene Name:	PEST proteolytic signal containing nuclear protein
Database Link:	NP_065090 Entrez Gene 57092 Human Q8WW12
Background:	PCNP may be involved in cell cycle regulation.
Synonyms:	DKFZp781I24156
Note:	Immunogen Sequence Homology: Dog: 100%; Pig: 100%; Rat: 100%; Horse: 100%; Human: 100%; Mouse: 100%; Bovine: 100%; Rabbit: 100%; Guinea pig: 100%



[View online »](#)

Product images:



WB Suggested Anti-PCNP Antibody Titration: 0.2-1 ug/ml; Positive Control: Hela cell lysate