

## Product datasheet for **TA344796**

### ANKMY2 Rabbit Polyclonal Antibody


#### Product data:

Product Type:	Primary Antibodies
Applications:	WB
Recommended Dilution:	WB
Reactivity:	Human
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	The immunogen for anti-ANKMY2 antibody: synthetic peptide directed towards the N terminal of human ANKMY2. Synthetic peptide located within the following region: MVHIKKGELTQEEKELLELVIGKGTVQEAGTLLSSKNVRVNCLDENGMTPL
Formulation:	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose. <i>Note that this product is shipped as lyophilized powder to China customers.</i>
Purification:	Affinity Purified
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Predicted Protein Size:	49 kDa
Gene Name:	ankyrin repeat and MYND domain containing 2
Database Link:	<a href="#">NP_064715</a> <a href="#">Entrez Gene 57037 Human</a> <a href="#">Q8IV38</a>
Background:	The function of this protein remains unknown.
Synonyms:	ZMYND20
Note:	Immunogen Sequence Homology: Pig: 100%; Rat: 100%; Horse: 100%; Human: 100%; Mouse: 100%; Bovine: 100%; Dog: 92%; Yeast: 90%; Guinea pig: 86%



[View online »](#)

## Product images:



WB Suggested Anti-ANKMY2 Antibody Titration:  
0.2-1 ug/ml; ELISA Titer: 1:312500; Positive  
Control: Human Muscle