

Product datasheet for **TA344868**


RSBN1 Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	WB
Recommended Dilution:	WB
Reactivity:	Human
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	The immunogen for anti-RSBN1 antibody: synthetic peptide directed towards the N terminal of human RSBN1. Synthetic peptide located within the following region: GGAVGPFKCVFVGEMAAQVGAVRVVRAVAAQEEDKEGKEKPHAGVSPRG
Formulation:	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Purification:	Affinity Purified
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Predicted Protein Size:	90 kDa
Gene Name:	round spermatid basic protein 1
Database Link:	NP_060834 Entrez Gene 54665 Human Q5VWQ0
Background:	The function of this protein remains unknown.
Synonyms:	ROSBIN
Note:	Immunogen Sequence Homology: Dog: 100%; Pig: 100%; Rat: 100%; Human: 100%; Mouse: 100%; Rabbit: 100%; Guinea pig: 100%; Bovine: 86%



[View online »](#)

Product images:

WB Suggested Anti-RSBN1 Antibody Titration:
0.2-1 ug/ml; ELISA Titer: 1:1562500; Positive
Control: Human Lung