

Product datasheet for **TA344901**

SMU1 Rabbit Polyclonal Antibody


Product data:

Product Type:	Primary Antibodies
Applications:	WB
Recommended Dilution:	WB
Reactivity:	Human
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	The immunogen for anti-SMU1 antibody: synthetic peptide directed towards the N terminal of human SMU1. Synthetic peptide located within the following region: KQTQPERYIHLENLLARSYFDPREAYPDGSSKEKRRAAIAQALAGEVSV
Formulation:	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Purification:	Affinity Purified
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Predicted Protein Size:	57 kDa
Gene Name:	DNA replication regulator and spliceosomal factor
Database Link:	NP_060695 Entrez Gene 55234 Human Q2TAY7
Background:	The function of this protein remains unknown.
Synonyms:	BWD; fSAP57; SMU-1
Note:	Immunogen Sequence Homology: Dog: 100%; Pig: 100%; Rat: 100%; Horse: 100%; Human: 100%; Mouse: 100%; Bovine: 100%; Rabbit: 100%; Zebrafish: 100%; Guinea pig: 100%



[View online »](#)

Product images:



WB Suggested Anti-SMU1 Antibody Titration: 0.2-1 ug/ml; Positive Control: HepG2 cell lysate