

Product datasheet for **TA345100**

FAM135B Rabbit Polyclonal Antibody


Product data:

Product Type:	Primary Antibodies
Applications:	WB
Recommended Dilution:	WB
Reactivity:	Human
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	The immunogen for anti-FAM135B antibody: synthetic peptide directed towards the middle region of human FAM135B. Synthetic peptide located within the following region: TLVSTGLWLMQKLKSGSLLQLTFRDNADLRKCFLYQLSQKTGLQYFKNV
Formulation:	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Purification:	Affinity Purified
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Predicted Protein Size:	156 kDa
Gene Name:	family with sequence similarity 135 member B
Database Link:	NP_056996 Entrez Gene 51059 Human Q49AJ0
Background:	The function of this protein remains unknown.
Synonyms:	C8ORFK32
Note:	Immunogen Sequence Homology: Dog: 100%; Pig: 100%; Rat: 100%; Horse: 100%; Human: 100%; Mouse: 100%; Bovine: 100%; Rabbit: 100%; Guinea pig: 100%



[View online »](#)

Product images:



WB Suggested Anti-FAM135B Antibody Titration:
0.2-1 ug/ml; ELISA Titer: 1:12500; Positive Control:
Human brain