

Product datasheet for **TA345269**

ZNF307 (ZKSCAN4) Rabbit Polyclonal Antibody

Product data:


Product Type:	Primary Antibodies
Applications:	WB
Recommended Dilution:	WB
Reactivity:	Human
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	The immunogen for anti-ZNF307 antibody: synthetic peptide directed towards the C terminal of human ZNF307. Synthetic peptide located within the following region: FTRNRS�IEHQKIHTGEKPYQCDCGKGFTRTSYLQHQQRSHVGGKKTLSQ
Formulation:	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose. <i>Note that this product is shipped as lyophilized powder to China customers.</i>
Purification:	Protein A purified
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Predicted Protein Size:	62 kDa
Gene Name:	zinc finger with KRAB and SCAN domains 4
Database Link:	NP_061983 Entrez Gene 387032 Human Q969J2
Background:	ZNF307 contains 1 SCAN box domain, 1 KRAB domain and 7 C2H2-type zinc fingers. It belongs to the Krueppel C2H2-type zinc-finger protein family and may be involved in transcriptional regulation.
Synonyms:	P1P373C6; ZNF307; ZNF427; ZSCAN36
Note:	Immunogen Sequence Homology: Rat: 100%; Horse: 100%; Human: 100%; Pig: 92%; Bovine: 92%; Dog: 86%; Mouse: 86%; Rabbit: 79%



[View online »](#)

Protein Families: Transcription Factors

Product images:



WB Suggested Anti-ZNF307 Antibody Titration:
1.25 ug/ml; ELISA Titer: 1: 312500; Positive
Control: Transfected 293T