


Product datasheet for **TA345305**

ZNF530 Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	WB
Recommended Dilution:	WB
Reactivity:	Human
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	The immunogen for anti-ZNF530 antibody: synthetic peptide directed towards the middle region of human ZNF530. Synthetic peptide located within the following region: CRECGKSFTTRKNHLIQHKT VHTGERPYECSECGKSFSQSSGLLRHRRVHV
Formulation:	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Purification:	Affinity Purified
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Predicted Protein Size:	69 kDa
Gene Name:	zinc finger protein 530
Database Link:	NP_065931 Entrez Gene 348327 Human Q6P9A1
Background:	ZNF530 belongs to the krueppel C2H2-type zinc-finger protein family. It contains 13 C2H2-type zinc fingers and 1 KRAB domain. ZNF530 may be involved in transcriptional regulation.
Synonyms:	KIAA1508
Note:	Immunogen Sequence Homology: Horse: 100%; Human: 100%; Dog: 86%; Pig: 86%; Rat: 86%; Bovine: 86%; Rabbit: 86%; Zebrafish: 86%; Goat: 79%; Sheep: 79%
Protein Families:	Transcription Factors


[View online »](#)

Product images:


WB Suggested Anti-ZNF530 Antibody Titration:
0.2-1 ug/ml; ELISA Titer: 1: 1562500; Positive
Control: HT1080 cell lysate