

Product datasheet for **TA345314**

Deformed Epidermal Autoregulatory Factor 1 (DEAF1) Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	WB
Recommended Dilution:	WB
Reactivity:	Human
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	The immunogen for anti-DEAF1 antibody: synthetic peptide directed towards the N terminal of human DEAF1. Synthetic peptide located within the following region: EEPVLSRDEDESEEDADSEAERETPRVTAVAVMAAEPGHMDMGAEALPGPD
Formulation:	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose. <i>Note that this product is shipped as lyophilized powder to China customers.</i>
Purification:	Affinity Purified
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Predicted Protein Size:	59 kDa
Gene Name:	DEAF1, transcription factor
Database Link:	NP_066288 Entrez Gene 10522 Human O75398
Background:	DEAF1 down-regulates transcription of those genes by binding to sequence with multiple copies of TTC[CG]G present in their own promoter and that of the HNRPA2B1 gene. DEAF1 binds to the retinoic acid response element (RARE) AGGGTTCACCGAAAGTTCA. DEAF1 activates the proenkephalin gene independently of promoter binding, probably through protein-protein interaction. When secreted, behaves as an inhibitor of cell proliferation, by arresting cells in the G0 or G1 phase.




[View online »](#)

Synonyms: MRD24; NUDR; SPN; ZMYND5

Note: Immunogen Sequence Homology: Pig: 100%; Rat: 100%; Human: 100%; Mouse: 100%; Yeast: 100%; Guinea pig: 100%; Dog: 93%; Bovine: 83%

Protein Families: Secreted Protein, Transcription Factors, Transmembrane

Product images:



WB Suggested Anti-DEAF1 Antibody Titration: 0.2-1 ug/ml; ELISA Titer: 1: 12500; Positive Control: Human brain

DEAF1



Rabbit Anti-DEAF1
Sample Type: Human Fetal Lung
Antibody Concentration: 1ug/mL

Hum. Fetal Lung

DEAF1



DEAF1

