

Product datasheet for **TA345331**

ZNF71 Rabbit Polyclonal Antibody

Product data:


Product Type:	Primary Antibodies
Applications:	IHC, WB
Recommended Dilution:	WB, IHC
Reactivity:	Human
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	The immunogen for anti-ZNF71 antibody: synthetic peptide directed towards the middle region of human ZNF71. Synthetic peptide located within the following region: RCGQCGKSFKNSSLTVHQRIHTGKPYRCGECGKTF SRNTNLTRHLRIH
Formulation:	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose. <i>Note that this product is shipped as lyophilized powder to China customers.</i>
Purification:	Affinity Purified
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Predicted Protein Size:	54 kDa
Gene Name:	zinc finger protein 71
Database Link:	NP_067039 Entrez Gene 58491 Human Q9NQZ8
Background:	ZNF71 belongs to the krueppel C2H2-type zinc-finger protein family and may be involved in transcriptional regulation.
Synonyms:	EZFIT
Note:	Immunogen Sequence Homology: Dog: 100%; Pig: 100%; Horse: 100%; Human: 100%; Rat: 86%; Mouse: 86%; Bovine: 86%; Rabbit: 86%; Guinea pig: 86%




[View online »](#)

Protein Families: Druggable Genome, Transcription Factors


Product images:




WB Suggested Anti-ZNF71 Antibody Titration: 1 ug/ml; Positive Control: PANC1 cell lysate




Anti-ZNF71 antibody IHC staining of human uterus. Immunohistochemistry of formalin-fixed, paraffin-embedded tissue after heat-induced antigen retrieval. Antibody concentration 5 ug/ml.



Anti-ZNF71 antibody IHC staining of human spleen. Immunohistochemistry of formalin-fixed, paraffin-embedded tissue after heat-induced antigen retrieval. Antibody concentration 5 ug/ml.



Anti-ZNF71 antibody IHC staining of human placenta. Immunohistochemistry of formalin-fixed, paraffin-embedded tissue after heat-induced antigen retrieval. Antibody concentration 5 ug/ml.



Anti-ZNF71 antibody IHC staining of human placenta. Immunohistochemistry of formalin-fixed, paraffin-embedded tissue after heat-induced antigen retrieval. Antibody concentration 5 ug/ml.