


Product datasheet for TA345336

PRDM13 Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	WB
Recommended Dilution:	WB
Reactivity:	Human
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	The immunogen for anti-PRDM13 antibody: synthetic peptide directed towards the N terminal of human PRDM13. Synthetic peptide located within the following region: HGAARAPATSVSADCCIPAGLRLGPVPGTFKLGKYLSDRREPGPKKKVRM
Formulation:	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Purification:	Affinity Purified
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Predicted Protein Size:	74 kDa
Gene Name:	PR domain 13
Database Link:	NP_067633 Entrez Gene 59336 Human Q9H4Q3
Background:	PRDM13 may be involved in transcriptional regulation.
Synonyms:	MU-MB-20.220; PFM10
Note:	Immunogen Sequence Homology: Dog: 100%; Human: 100%; Bovine: 100%; Pig: 93%; Rat: 93%; Horse: 93%; Guinea pig: 86%; Rabbit: 79%


[View online »](#)

Product images:

WB Suggested Anti-PRDM13 Antibody Titration:
0.2-1 ug/ml; ELISA Titer: 1: 312500; Positive
Control: Human heart