

Product datasheet for **TA345343**

ZNF70 Rabbit Polyclonal Antibody


Product data:

Product Type:	Primary Antibodies
Applications:	IHC, WB
Recommended Dilution:	WB, IHC
Reactivity:	Human
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	The immunogen for anti-ZNF70 antibody: synthetic peptide directed towards the C terminal of human ZNF70. Synthetic peptide located within the following region: GKAFRHRSALIEHYKTHTREKPYVCNLCGKSFRGSSHLIRHQIHSGEKL
Formulation:	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Purification:	Protein A purified
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Predicted Protein Size:	51 kDa
Gene Name:	zinc finger protein 70
Database Link:	NP_068735 Entrez Gene 7621 Human Q9UC06
Background:	The function of the ZNF70 has not yet been determined
Synonyms:	Cos17
Note:	Immunogen Sequence Homology: Dog: 100%; Rat: 100%; Human: 100%; Pig: 93%; Horse: 93%; Bovine: 86%; Rabbit: 83%; Zebrafish: 83%; Mouse: 79%; Guinea pig: 77%
Protein Families:	Transcription Factors




[View online »](#)

Product images:



WB Suggested Anti-ZNF70 Antibody Titration: 5.0 ug/ml; ELISA Titer: 1: 312500; Positive Control: HepG2 cell lysate



Human Stomach