

Product datasheet for **TA345377**

HOXB4 Rabbit Polyclonal Antibody


Product data:

Product Type:	Primary Antibodies
Applications:	WB
Recommended Dilution:	WB
Reactivity:	Human
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	The immunogen for anti-HOXB4 antibody: synthetic peptide directed towards the N terminal of human HOXB4. Synthetic peptide located within the following region: FLINSNYVDPKFPPEEYSQSDYLPSDHSPGYAGGQRRESSFQPEAGFG
Formulation:	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose. <i>Note that this product is shipped as lyophilized powder to China customers.</i>
Purification:	Protein A purified
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Predicted Protein Size:	28 kDa
Gene Name:	homeobox B4
Database Link:	NP_076920 Entrez Gene 3214 Human P17483



[View online »](#)

Background:	HOXB4 is a nuclear protein with a homeobox DNA-binding domain. The protein functions as a sequence-specific transcription factor that is involved in development. Intracellular or ectopic expression of this protein expands hematopoietic stem and progenitor cells in vivo and in vitro, making it a potential candidate for therapeutic stem cell expansion. This gene is a member of the Antp homeobox family and encodes a nuclear protein with a homeobox DNA-binding domain. It is included in a cluster of homeobox B genes located on chromosome 17. The encoded protein functions as a sequence-specific transcription factor that is involved in development. Intracellular or ectopic expression of this protein expands hematopoietic stem and progenitor cells in vivo and in vitro, making it a potential candidate for therapeutic stem cell expansion.
Synonyms:	HOX-2.6; HOX2; HOX2F
Note:	Immunogen Sequence Homology: Dog: 100%; Pig: 100%; Rat: 100%; Human: 100%; Mouse: 100%; Bovine: 100%
Protein Families:	ES Cell Differentiation/IPS, Stem cell relevant signaling - DSL/Notch pathway, Transcription Factors

Product images:

WB Suggested Anti-HOXB4 Antibody Titration: 1.0 ug/ml; ELISA Titer: 1: 62500; Positive Control: Jurkat cell lysate