

## Product datasheet for **TA345385**

### ZSCAN5 (ZSCAN5A) Rabbit Polyclonal Antibody


#### Product data:

Product Type:	Primary Antibodies
Applications:	WB
Recommended Dilution:	WB
Reactivity:	Human
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	The immunogen for Anti-ZSCAN5A antibody is: synthetic peptide directed towards the C-terminal region of Human ZSCAN5A. Synthetic peptide located within the following region: KVFTYRGSLKEHQRIHSGEKPYKCSKCPRAFSRLKLLRRHQKTHPEATSQ
Formulation:	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Purification:	Affinity Purified
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Predicted Protein Size:	41 kDa
Gene Name:	zinc finger and SCAN domain containing 5A
Database Link:	<a href="#">NP_077279</a> <a href="#">Entrez Gene 79149 Human</a> <a href="#">Q9BUG6</a>
Background:	The function of this protein remains unknown.
Synonyms:	ZNF495; ZSCAN5
Note:	Immunogen Sequence Homology: Human: 100%
Protein Families:	Transcription Factors



[View online »](#)

## Product images:



Host: Rabbit; Target Name: ZSCAN5A; Sample Tissue: Liver tumor lysates; Antibody Dilution: 1.0 ug/ml